

1	1	10	20	BM-HABP Fig. 4
				TSG-6.PRO
21	3	30	40	BM-HABP Fig. 4
				TSG-6.PRO
41	8	50	60	BM-HABP Fig. 4
				TSG-6.PRO
49	24	70	80	BM-HABP Fig. 4
				TSG-6.PRO
64	44	90	100	BM-HABP Fig. 4
				TSG-6.PRO
84	64	110	120	BM-HABP Fig. 4
				TSG-6.PRO
104	84	130	140	BM-HABP Fig. 4
				TSG-6.PRO
124	104	150	160	BM-HABP Fig. 4
				TSG-6.PRO
144	124	170	180	BM-HABP Fig. 4
				TSG-6.PRO
159	144	190	200	BM-HABP Fig. 4
				TSG-6.PRO
174	164	210	220	BM-HABP Fig. 4
				TSG-6.PRO
182	164	230	240	BM-HABP Fig. 4
				TSG-6.PRO
197	204	250	260	BM-HABP Fig. 4
				TSG-6.PRO
217	206	270	280	BM-HABP Fig. 4
				TSG-6.PRO
235	226	290	300	BM-HABP Fig. 4
				TSG-6.PRO
255	236	310	320	BM-HABP Fig. 4
				TSG-6.PRO
275	249	330	340	BM-HABP Fig. 4
				TSG-6.PRO
295	263	350	360	BM-HABP Fig. 4
				TSG-6.PRO
315	270	370	380	BM-HABP Fig. 4
				TSG-6.PRO
335	270	390		BM-HABP Fig. 4
				TSG-6.PRO

Percent Identity			
Divergence		1	2
	1		31.6
	2	100.0	
		1	2
BM-HABP Fig.4.PRO			
TSG-6.PRO			



Blast 2 Sequences results


[PubMed](#)[Entrez](#)[BLAST](#)[OMIM](#)[Taxonomy](#)[Structure](#)

BLAST 2 SEQUENCES RESULTS VERSION BLASTP 2.1.2 [Oct-19-2000]

Matrix: **BLOSUM62** gap open: **11** gap extension: **1**
x_dropoff: **50** expect: **10.000** wordsize: **3** Filter ☒ **Align**

Sequence 1 lc1lsèq_1 Length 353 (1 .. 353)

Sequence 2 lc1lsq_2 Length 275 (1 .. 275)



NOTE: The statistics (bitscore and expect value) is calculated based on the size of nr database

Score = 107 bits (264), Expect = 5e-21
Identities = 45/104 (43%), Positives = 61/104 (58%)

Query: 52 DTTVGVFHLRSPLGQYKLTFDKAREACANEAATMATYNQLSYXQKAKYHLCSAGWLETGR 111
+ GV+H + G+YKLT+ +A+ C E +ATY QL +K +H+C+AGW+ GR
Sbjct: 32 EQAAGVYHREARAGRYKLTAEAKAVCEFEGRRLATYKQLEAARKIGFHVCAAGWMAKGR 91

Query: 112 VAYPTAFASQNCGSGVVGIVDYGPRPNKSEMWDVFCYRMKDVCNC 155
V YP NCG G GI+DYG R N+SE WD +CY C
Sbjct: 92 VGYPPIVKPGPNCGFGKTGIIDYGIRLNRSERWDAYCYNPHAKEC 135

CPU time: 0.15 user secs. 0.03 sys. secs 0.18 total secs.

Gapped

Lambda	K	H
0.321	0.138	0.429

Gapped

Lambda	K	H
0.270	0.0470	0.230

Matrix: BLOSUM62

Gap Penalties: Existence: 11, Extension: 1

Number of Hits to DB: 748

Number of Sequences: 0

Number of extensions: 55

Number of successful extensions: 1

Number of sequences better than 10.0: 1

Number of HSP's better than 10.0 without gapping: 1

Number of HSP's successfully gapped in prelim test: 0

Number of HSP's that attempted gapping in prelim test: 0

Number of HSP's gapped (non-prelim): 1

length of query: 353
length of database: 3,171,650,076
effective HSP length: 58
effective length of query: 295
effective length of database: 3,171,650,018
effective search space: 935636755310
effective search space used: 935636755310
T: 9
A: 40
X1: 16 (7.4 bits)
X2: 128 (49.9 bits)
X3: 128 (49.9 bits)
S1: 41 (21.9 bits)
S2: 83 (36.7 bits)

EXHIBIT C

56 GVFHLRSPLGQYKLTFDKAREACANEAATMATYNQLSYXQKAKYHLC SAGWLETGRVAYP **-SEQ ID NO:11**
|||||
1063 GVFHLRSPLGQYKLTFDKAREACANEAATMATYNQLSYAQKAKYHLC SAGWLETGRVAYP **-human HARE**

115 TAFASQNCGSGVVGIVDYGPRPNKSEMWDVFCYR 149
|||||
1122 TAFASQNCGSGVVGIVDYGPRPNKSEMWDVFCYR 1156

CURRENT PROTOCOLS IN MOLECULAR BIOLOGY

VOLUME 1

EDITORIAL BOARD

Frederick M. Ausubel
Massachusetts General Hospital & Harvard Medical School

Roger Brent
Massachusetts General Hospital & Harvard Medical School

Robert E. Kingston
Massachusetts General Hospital & Harvard Medical School

David D. Moore
Massachusetts General Hospital & Harvard Medical School

J.G. Seidman
Harvard Medical School

John A. Smith
University of Alabama

Kevin Struhl
Harvard Medical School

GUEST EDITORS

Lisa M. Albright
DNA Sequencing

Donald M. Coen
Harvard Medical School
Polymerase Chain Reaction

Ajit Varki
University of California San Diego
Glycoproteins

SERIES EDITOR

Virginia Benson Chanda



John Wiley & Sons, Inc.

CORE 13 (S36)

Copyright © 1994–1997 by John Wiley & Sons, Inc.

Copyright © 1987–1994 by Current Protocols

All rights reserved. Published simultaneously in Canada.

Reproduction or translation of any part of this work beyond that permitted by Section 107 or 108 of the 1976 United States Copyright Act without the permission of the copyright owner is unlawful. Requests for permission or further information should be addressed to the Permissions Department, John Wiley & Sons, Inc.

While the authors, editors, and publisher believe that the specification and usage of reagents, equipment, and devices, as set forth in this book, are in accord with current recommendations and practice at the time of publication, they accept no legal responsibility for any errors or omissions, and make no warranty, express or implied, with respect to material contained herein. In view of ongoing research, equipment modifications, changes in governmental regulations, and the constant flow of information relating to the use of experimental reagents, equipment, and devices, the reader is urged to review and evaluate the information provided in the package insert or instructions for each chemical, piece of equipment, reagent, or device for, among other things, any changes in the instructions or indication of usage and for added warnings and precautions. This is particularly important in regard to new or infrequently employed chemicals or experimental reagents.

Library of Congress Cataloging in Publication Data:

Current protocols in molecular biology. 3 vols.

1. Molecular biology—Technique. 2. Molecular biology—Laboratory manuals. I. Ausbel, Frederick M.

QH506.C87 1987 574.873028 87-21033

ISBN 0-471-50338-X

Printed in the United States of America

20 19 18 17 16 15 14 13

- 1c. *Harsh treatment:* Pour several hundred milliliters of boiling 0.1% SDS onto the membrane. Cool to room temperature.

If a membrane is to be reprobed, it must not be allowed to dry out between hybridization and stripping. If it becomes dry, the probe may bind to the matrix.

2. Place membrane on a sheet of dry Whatman 3MM filter paper and blot excess liquid with a second sheet. Wrap the membrane in plastic wrap and set up an autoradiograph.

If signal is still seen after autoradiography, rewash using harsher conditions.

3. The membrane can now be rehybridized. Alternatively, it can be dried and stored for later use.

Membranes can be stored dry between Whatman 3MM paper for several months at room temperature. For long-term storage, place the membranes in a desiccator at room temperature or 4°C.

REAGENTS AND SOLUTIONS

Aqueous prehybridization/hybridization (APH) solution

5× SSC (APPENDIX 2)

5× Denhardt solution (APPENDIX 2)

1% (w/v) SDS

Add 100 µg/ml denatured salmon sperm DNA (see below) just before use

Alternatives to Denhardt solution and denatured salmon sperm DNA as blocking agents are listed in Table 2.10.5 (see discussion in critical parameters).

Denatured salmon sperm DNA

Dissolve 10 mg Sigma type III salmon sperm DNA (sodium salt) in 1 ml water. Pass vigorously through a 17-G needle 20 times to shear the DNA. Place in a boiling water bath for 10 min, then chill. Use immediately or store at -20°C in small aliquots. If stored, reheat to 100°C for 5 min and chill on ice immediately before using.

Formamide prehybridization/hybridization (FPH) solution

5× SSC (APPENDIX 2)

5× Denhardt solution (APPENDIX 2)

50% (w/v) formamide

1% (w/v) SDS

Add 100 µg/ml denatured salmon sperm DNA (see above) just before use

Alternatives to Denhardt solution and denatured salmon sperm DNA as blocking agents are listed in Table 2.10.5 (see discussion in critical parameters).

Commercial formamide is usually satisfactory for use. If the liquid has a yellow color, deionize as follows: add 5 g of mixed-bed ion-exchange resin [e.g., Bio-Rad AG 501-X8 or 501-X8(D) resins] per 100 ml formamide, stir at room temperature for 1 hr, and filter through Whatman no. 1 paper.

CAUTION: Formamide is a teratogen. Handle with care.

Labeling buffer

200 mM Tris-Cl, pH 7.5

30 mM MgCl₂

10 mM spermidine

Mild stripping solution

5 mM Tris-Cl, pH 8.0

2 mM EDTA

0.1× Denhardt solution (APPENDIX 2)

SDS electrophoresis buffer, 5×

15.1 g Tris base

72.0 g glycine

5.0 g SDS

H₂O to 1000 ml

Dilute to 1× or 2× for working solution, as appropriate

Do not adjust the pH of the stock solution, as the solution is pH 8.3 when diluted. Store at 0° to 4° C until use (up to 1 month).

SED (standard enzyme diluent)

20 mM Tris-Cl, pH 7.5

500 µg/ml bovine serum albumin (Pentax Fraction V)

10 mM 2-mercaptoethanol

Store up to 1 month at 4° C

Sodium acetate, 3 M

Dissolve 408 g sodium acetate-3H₂O in 800 ml H₂O

Add H₂O to 1 liter

Adjust pH to 4.8 or 5.2 (as desired) with 3 M acetic acid

Sodium acetate buffer, 0.1 M

Solution A: 11.55 ml glacial acetic acid/liter (0.2 M).

Solution B: 27.2 g sodium acetate (NaC₂H₃O₂·3H₂O)/liter (0.2 M).

Referring to Table A.2.2 for desired pH, mix the indicated volumes of solutions A and B, then dilute with H₂O to 100 ml. (See Potassium acetate buffer recipe for further details.)

Sodium phosphate buffer, 0.1 M

Solution A: 27.6 g NaH₂PO₄·H₂O per liter (0.2 M).

Solution B: 53.65 g Na₂HPO₄·7H₂O per liter (0.2 M).

Referring to Table A.2.3 for desired pH, mix the indicated volumes of solutions A and B, then dilute with H₂O to 200 ml. (See Potassium phosphate buffer recipe for further details.)

SSC (sodium chloride/sodium citrate), 20×

3 M NaCl (175 g/liter)

0.3 M Na₃citrate-2H₂O (88 g/liter)

Adjust pH to 7.0 with 1 M HCl

STE buffer

10 mM Tris-Cl, pH 7.5

10 mM NaCl

1 mM EDTA, pH 8.0

TAE (Tris/acetate/EDTA) electrophoresis buffer

50× stock solution:

242 g Tris base

57.1 ml glacial acetic acid

37.2 g Na₂EDTA-2H₂O

H₂O to 1 liter

Working solution, pH ~8.5:

40 mM Tris-acetate

2 mM Na₂EDTA-2H₂O

TBE (Tris/borate/EDTA) electrophoresis buffer

10× stock solution, 1 liter:

108 g Tris base (890 mM)

55 g boric acid (890 mM)

40 ml 0.5 M EDTA, pH 8.0 (20 mM)